

Shenzhen Pango Medical Electronics Co., Ltd.
Main Site: Building 2, No. 25 Fenghuang Road, Industrial Zone, Xixiang First Village, Henggang Street, Longgang District, Shenzhen, 518115 Guangdong, P.R. China.
Additional site 1: 2-4 Floor, No.5 Shanzhang Rd., Xixiang Village Henggang Town, Longgang District, Shenzhen, 518115 Guangdong, P.R. China.
Tel:+86-755-33825988 Fax:+86-755-33825989

Lotus NL B.V.
Address: Konings Julianenplein 10, 1e Verd, 2555AA, The Hague, Netherlands.
Tel: +31644168999

Date: 2023-04-13 Rev:A/2

ELECTRONIC BLOOD PRESSURE MONITOR

Instruction Manual

MODEL: PG-800A52

TABLE OF CONTENTS	
INTRODUCTION.....	2
NOTES ON SAFETY.....	2
ABOUT BLOOD PRESSURE.....	6
PRECAUTIONS BEFORE USE.....	8
FEATURES OF THE PRODUCT.....	9
PARTS IDENTIFICATION.....	10
INSERT OR REPLACE BATTERIES.....	11
TIME AND VOICE ON/OFF OF SYSTEM SETUP.....	12
UNIT CONVERSION mmHg/kPa DISPLAY.....	13
WHO BLOOD PRESSURE CLASSIFICATION DISPLAY.....	14
ATTACHING THE WRIST CUFF.....	15
HOW TO MEASURE BLOOD PRESSURE.....	16
CLEAN AND MAINTENANCE.....	17
SPECIFICATIONS.....	19
TROUBLESHOOTING.....	20
CALIBRATION METHOD.....	29

INTRODUCTION	
<p>The Monitor uses the oscillometric method of blood pressure measurement. Intended for use by medical professionals or at home to monitor and display diastolic, systolic blood pressure and pulse rate on adult.</p> <p>With an air wrist cuff buckled around one's wrist according to the instructions in the " ATTACHING THE WRIST CUFF".</p> <p>The expected life of the product is 5 years.</p> <p>The product complies with the electromagnetic compatibility requirement of IEC 60601-1-2 and safety standards of IEC 60601-1 and performance of IEC 80601-2-30 as specified in Regulation (EU)2017/745.</p>	
<p>* The warning signs and samples icons shown here are listed for your safe and correct use of the unit, so as to prevent injuries or damages to the device.</p> <p>* The icons and meanings are as follow.</p> <p>Examples of signs</p> <p>The ☹ icon indicates prohibitions (what you should not do). Matters involving actual prohibitions are indicated by text or pictures in or near ☹. The left icon refers to "general prohibition".</p>	

NOTES ON SAFETY	
<p>The ☹ icon indicates something that is compulsory (what must always be observed).</p> <p>Matters involving actual compulsory actions are indicated by text or pictures in or near ☹. The left icon refers to "general prohibition".</p> <p>The ☹ icon indicates something can't be disassembled or "Don't disassemble" Matters involving actual compulsory actions are indicated by text or pictures in or near ☹. The left icon refers to "general prohibition".</p>	
<p>Type BF Applied part</p> <p>IP Classification: IP22</p> <p>Caution</p> <p>Please refer to the instructions for use</p> <p>Consult instruction for use</p> <p>Indicates a medical device that needs to be protected from moisture.</p> <p>Contact its local authorities to determine the proper method of disposal of potentially bio hazardous parts and accessories.</p>	<p>The following symbol indicates that the device is MR-unsafe:</p>
<p>Do not measure your blood pressure over 6 times each day.</p> <p>Do not apply the cuff over a wound as this can cause further injury.</p> <p>Do not measure on the wrist which is on the side of a mastectomy, otherwise it could cause injury.</p> <p>Observe the air pressure value from the LCD display.</p> <p>When measuring, it could not exceed 280 mmHg, otherwise Please press "on/off" button to stop</p> <p>Do not use force to bend the wrist cuff or the air tube.</p> <p>Do not knock or drop the main unit.</p> <p>Always use the specified accessories in the manual, the use of other parts not approved by the manufacturer may cause faults or injuries.</p> <p>For service information, parts list etc., please contact the dealer.</p>	

ABOUT BLOOD PRESSURE	
<p>1. What is blood pressure?</p> <p>Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts. Diastolic pressure occurs when the heart expands.</p> <p>Blood pressure is measured in millimeters of mercury (mmHg). One's natural blood pressure is represented by the fundamental pressure, which is measured first thing in the morning while one is still at rest and before eating.</p> <p>2. What is hypertension and how is it controlled?</p> <p>Hypertension, an abnormally high arterial blood pressure, if left unattended, can cause many health problems including stroke and heart attack. Hypertension can be controlled by altering lifestyle, avoiding stress and with medication under a doctor's supervision. To prevent hypertension or keep it under control:</p>	<p>The PATIENT is an intended OPERATOR.</p> <p>-Not servicing and maintenance while the ME EQUIPMENT is in use.</p> <p>-The user can maintain the product, the maintenance method is described in the maintenance instructions of manual.</p> <p>-Stop using the equipment immediately, if it is in contact with water.</p>
<p>Make sure there is no connection tubing kinking before start measuring to avoid any injury to patient.</p> <p>For any patient, do not measure more than 3 times continuously, it should be at least above 5 minutes of interval rest between any two measurements, otherwise will cause extravasated blood.</p> <p>Do not measure your blood pressure over 6 times each day.</p> <p>Do not apply the cuff over a wound as this can cause further injury.</p> <p>Do not measure on the wrist which is on the side of a mastectomy, otherwise it could cause injury.</p> <p>Observe the air pressure value from the LCD display.</p> <p>When measuring, it could not exceed 280 mmHg, otherwise Please press "on/off" button to stop</p> <p>Do not use force to bend the wrist cuff or the air tube.</p> <p>Do not knock or drop the main unit.</p> <p>Always use the specified accessories in the manual, the use of other parts not approved by the manufacturer may cause faults or injuries.</p> <p>For service information, parts list etc., please contact the dealer.</p>	

PRECAUTIONS BEFORE USE	
<p>1. If you are taking medication, consult with your doctor to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your doctor.</p>	<p>hypertensive individuals, variations are even more pronounced.</p> <p>Normally, the blood pressure rises while at work or play and falls to its lowest levels during sleep. So, do not be overly concerned by the results of one measurement.</p> <p>Take measurements at the same time every day using the procedure described in this manual, and know your normal blood pressure. Many readings give a more comprehensive blood pressure history.</p> <p>Be sure to note date and time when recording your blood pressure. Consult your doctor to interpret your blood pressure data.</p>

FEATURES OF THE PRODUCT	
<p>1. Memory can store 90 measurements.</p> <p>2. Large and clear LCD display.</p>	<p>2. For people with irregular or unstable peripheral circulation problems due to diabetes, liver disease, hardening of the arteries, etc., there may be fluctuation in blood pressure values measured at the upper arm versus at the wrist.</p> <p>3. Measurements may be impaired if this device is used near televisions, microwave ovens, X-ray, mobile phone equipment or other devices with strong electrical fields. To prevent such interference, use the monitor at a sufficient distance from such devices or turn them off.</p> <p>4. Before using, should wash your hands.</p> <p>5. Do not measure on the arm which simultaneously used monitoring ME Equipment, otherwise it could cause loss of function.</p> <p>6 Consult your doctor if the unexpected readings are obtained, also please refer to "Trouble shooting" of the manual.</p> <p>7. The reading is probably a little lower than measured in the hospital due to the steady mode at home.</p> <p>8. Cuff pressure range 0-299mmHg</p>

PARTS IDENTIFICATION	
<p>LCD Display: Date and Time, Systolic Blood pressure, Diastolic blood pressure, Pulse/beat/min</p> <p>SET Button</p> <p>ON/OFF Button</p> <p>Memory Button</p> <p>Battery Cover</p> <p>Wrist Cuff</p> <p>Accessory: Manual, CASE</p>	<p>SYMBOLS ON DISPLAY</p> <p>MEMO/DATE</p> <p>VALUE OF SYSTOLIC</p> <p>VALUE OF DIASTOLIC</p> <p>VALUE OF PULSE/BEAT</p> <p>MINI/HEART RATE</p> <p>MEASUREMENT DETECTION</p>

INSERT OR REPLACE BATTERIES	
<p>1. Remove the battery cover.</p> <p>2. Insert new batteries into the battery compartment as shown, taking care that the polarities(+) and (-)are correct.</p> <p>3. Close the battery cover. Use only LR03, AAA batteries.</p> <p>CAUTION</p> <p>• Insert the batteries as shown in the battery compartment. If not, the device will not work.</p> <p>• When LOW BATTERY mark) blinks in the display, replace all batteries with new ones. Do not mix old and new batteries. It may shorten the battery life, or cause the device to malfunction.</p> <p>• ☹ (LOW BATTERY mark) does not appear when the batteries run out.</p> <p>• Please ensure to distinguish positive polar "+" and negative polar "-" of batteries when replacing batteries.</p>	<p>• Batteries, which have fluid on surface or be modified, can not be inserted into the products.</p> <p>• Battery short circuit must be prevented.</p> <p>• Battery life varies with the ambient temperature and may be shorten at low temperatures.</p> <p>• The batteries may leak and cause a malfunction.</p> <p>• Use the specified batteries only. The batteries provided with the device are for testing monitor performance and may have a shorter life.</p> <p>• Used batteries may leak and damage the main unit. Pleases observe the following points.</p> <p>* If you are not going to use the unit for a long period of time (approximately three months or more), remove the batteries.</p> <p>* Replace worn batteries with their polarities in the correct direction.</p> <p>TIME AND VOICE ON/OFF OF SYSTEM SETUP</p> <p>1. Press "SET" key to Time display.</p> <p>2. In the off state, Press and hold "SET" key until the year number displays in flashes on LCD to enter setting mode.</p> <p>3. Press "MEM" key to adjust the year, then press "SET" key again to save your setting and enter the month setting mode.</p>

UNIT CONVERSION mmHg/kPa DISPLAY	
<p>The goods have mm Hg(mmHg), kPa (kPa) two kinds of blood pressure display units/(mmHg factory to express).</p> <p>Press "ON / OFF" button for 10 seconds to display unit switching interface, then press "MEM" key to select mmHg / KPa, press "ON / OFF" button to exit.</p>	<p>The units will be chosen by the above shows mmHg/kPa after decontrol, After the normal boot unit values are shown as blood pressure.</p> <p>Also select memory unit value changes.</p>

ATTACHING THE WRIST CUFF	
<p>1. Fastening the wrist cuff</p> <p>1) Wrap the wrist cuff around your wrist about (1-2)cm above your hand as shown in the figure at the right.</p> <p>2) Fasten the wrist cuff tightly by using the Velcro Strip.</p> <p>For proper measurements, fasten the wrist cuff tightly and measure on a bare wrist.</p> <p>2. How to take proper measurements</p> <p>For best accuracy in blood pressure measurement:</p> <ul style="list-style-type: none"> • Sit comfortably at a table. Rest your wrist on the table. • Relax for about 5 to 10 minutes before measurement. • Raise your hand so that the wrist cuff is at the same level as your heart. • Remain still and keep quiet during measurement. • Do not measure left after physical exercise or a bath. • Measure your blood pressure at about the same time every day. 	<p>3. WHO blood pressure classification display.</p> <p>4. Easy to use, Press a button to automatically measure, record the measurement values and flashes on LCD to enter setting mode.</p> <p>5. Automatically turns off (within 1 minute) to save power.</p>

HOW TO MEASURE BLOOD PRESSURE	
<p>1. Fasten the wrist cuff according to the instructions in "ATTACHING THE WRIST CUFF".</p> <p>2. The clean steps for the cuff is provided as following</p> <p>• Completely wipe the inner side (the side that contacts skin) of the cuff with a soft cloth lightly moistened with 75% Ethyl alcohol 3 times.</p> <p>• Replace the soft cloth after each wipe.</p> <p>• Then air dry the cuff.</p> <p>CAUTION</p> <p>• Do not submerge the device or any of the components in water.</p> <p>• Do not subject the monitor to extreme hot or cold temperatures, humidity or direct sunlight.</p> <p>• Store the device and the components in a clean, safe location.</p> <p>• Do not subject the monitor to strong shocks, such as dropping the unit on the floor.</p> <p>• Remove the batteries if the unit will not be used for three months or longer. Always replace all the batteries with new ones at the same time.</p> <p>• This product is designed for use over an extended period of time; however, it is generally recommended that it be inspected and calibrated every two years to ensure proper function and performance.</p> <p>• See the Calibration Method for more details.</p>	<p>4. Do not submerge the device or any components in water</p> <p>5. Store the device and the components in a clean and safe location.</p> <p>6. The clean steps for the cuff is provided as following</p> <p>• Completely wipe the inner side (the side that contacts skin) of the cuff with a soft cloth lightly moistened with 75% Ethyl alcohol 3 times.</p> <p>• Replace the soft cloth after each wipe.</p> <p>• Then air dry the cuff.</p>

CLEAN AND MAINTENANCE	
<p>1. Keep this device in the case provided with the device when you do not use it.</p> <p>2. Do not fold the arm cuff too tightly.</p> <p>3. Clean the monitor with a soft dry cloth. Do not use any cleaning solution.</p>	<p>18</p>

SPECIFICATIONS	
<p>Measuring Method: Oscillometric Measurement</p> <p>Indication: Digital LCD display</p> <p>Measuring Range: Pressure:(30~280)mmHg Pulse:(40~199)Beat/min</p> <p>Accuracy: Static Pressure: ±3mmHg Pulse: ±5%</p> <p>Memory: 90 Memories</p> <p>Power supply: 2x1.5V Batteries(LR03 or AAA) use alkaline battery measure above 200 times.</p> <p>Operating condition: +5℃~+40℃, 15%RH~93%RH Atmospheric pressure: 70kPa~106kPa</p> <p>Storage condition: -20℃~+55℃, 0%RH~93%RH Atmospheric pressure:50kPa~106kPa</p> <p>Dimensions: Approx: 72(W)X67(H)X30(D)mm</p> <p>Weight: Approx: 130g, excluding batteries</p> <p>Classification: Type BF</p> <p>Wrist circumference: (13.5~19.5)cm</p>	<p>1. Type of protection against electric shock: INTERNALLY POWERED EQUIPMENT.</p> <p>2. Degree of protection against electric shock: TYPE BF APPLIED PART.</p> <p>3. Mode of operation: CONTINUOUS OPERATION.</p> <p>4. Equipment not suitable for category AP&APG equipment use in presence.</p> <p>STATEMENT</p> <p>the system might not meet its performance specifications if stored or used outside the temperature and humidity as mentioned below.</p> <p>Operating conditions: +5℃~+40℃, 15%RH~93%RH 70kPa~106kPa</p> <p>Storage conditions: -20℃~+55℃, 0%RH~93%RH</p>

TROUBLE SHOOTING	
<p>If you have trouble in using the unit please check the following points first.</p>	<p>1. Error DISPLAY</p> <p>2. Possible CAUSE</p> <p>3. HOW TO CORRECT</p>
<p>Nothing is displayed</p> <p>When you push the POWER button or Battery icon flash</p>	<p>No battery installation</p> <p>Battery worn out</p> <p>The polarities of batteries placed wrongly</p>

Guidance and manufacturer's declaration – electromagnetic immunity	
<p>The Model PG-800A52 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A52 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.</p>	<p>Appendix 1 Guidance and Manufacturer Declaration Tables</p> <p>Guidance and manufacturer's declaration – electromagnetic emissions</p> <p>The Model PG-800A52 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A52 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.</p>
<p>Immunity test</p> <p>Electrostatic discharge (ESD) IEC 61000-4-2</p> <p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>Emissions</p> <p>RF emissions CISPR 11</p> <p>RF emissions CISPR 11</p> <p>Harmonic emissions IEC 61000-3-2</p> <p>Voltage fluctuations/flicker emissions IEC 61000-3-3</p>

Guidance and manufacturer's declaration – electromagnetic immunity	
<p>The Model PG-800A52 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A52 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.</p>	<p>Appendix 1 Guidance and Manufacturer Declaration Tables</p> <p>Guidance and manufacturer's declaration – electromagnetic immunity</p> <p>The Model PG-800A52 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A52 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.</p>
<p>Immunity test</p> <p>Electrostatic discharge (ESD) IEC 61000-4-2</p> <p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>Emissions</p> <p>RF emissions CISPR 11</p> <p>RF emissions CISPR 11</p> <p>Harmonic emissions IEC 61000-3-2</p> <p>Voltage fluctuations/flicker emissions IEC 61000-3-3</p>

Guidance and manufacturer's declaration – electromagnetic immunity	
<p>The Model PG-800A52 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A52 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.</p>	<p>Appendix 1 Guidance and Manufacturer Declaration Tables</p> <p>Guidance and manufacturer's declaration – electromagnetic immunity</p> <p>The Model PG-800A52 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A52 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.</p>
<p>Immunity test</p> <p>Electrostatic discharge (ESD) IEC 61000-4-2</p> <p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>Emissions</p> <p>RF emissions CISPR 11</p> <p>RF emissions CISPR 11</p> <p>Harmonic emissions IEC 61000-3-2</p> <p>Voltage fluctuations/flicker emissions IEC 61000-3-3</p>

Guidance and manufacturer's declaration – electromagnetic immunity	
<p>The Model PG-800A52 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A52 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.</p>	<p>Appendix 1 Guidance and Manufacturer Declaration Tables</p> <p>Guidance and manufacturer's declaration – electromagnetic immunity</p> <p>The Model PG-800A52 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A52 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.</p>
<p>Immunity test</p> <p>Electrostatic discharge (ESD) IEC 61000-4-2</p> <p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>Emissions</p> <p>RF emissions CISPR 11</p> <p>RF emissions CISPR 11</p> <p>Harmonic emissions IEC 61000-3-2</p> <p>Voltage fluctuations/flicker emissions IEC 61000-3-3</p>

Guidance and manufacturer's declaration – electromagnetic immunity	
<p>The Model PG-800A52 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A52 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.</p>	<p>Appendix 1 Guidance and Manufacturer Declaration Tables</p> <p>Guidance and manufacturer's declaration – electromagnetic immunity</p> <p>The Model PG-800A52 Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800A52 Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.</p>
<p>Immunity test</p> <p>Electrostatic discharge (ESD) IEC 61000-4-2</p> <p>Power frequency (50/60 Hz) magnetic field IEC 61000-4-8</p>	<p>Emissions</p> <p>RF emissions CISPR 11</p> <p>RF emissions CISPR 11</p> <p>Harmonic emissions IEC 61000-3-2</p> <p>Voltage fluctuations/flicker emissions IEC 61000-3-3</p>