

ELECTRONIC BLOOD PRESSURE MONITOR



Instruction Manual

MODEL: PG-800AD
PG-800AD-1



Shenzhen Pango Medical Electronics Co., Ltd
Main Site: Building 2, No. 25 Fenghuang Road, Industrial Zone, Xixiang First Village, Henggang Street, Longgang District, Shenzhen, 518115 Guangdong, P. R. China.
Additional site 1: 2-4 Floor, No. 5 Shanzhuang Rd., Xixiang Village, Henggang Town, Longgang District, Shenzhen, 518115 Guangdong, P. R. China.
Tel: +86-755-33825983 Fax: +86-755-33825989

Lotus N.L.B.V.
Address: Koningsn Julianenplein 10, 1e Veer, 2959AA, The Hague, Netherlands.
Tel: +31(0)441609999

Date: 2023-04-13 Rev:A/2

TABLE OF CONTENTS	
INTRODUCTION	2
NOTES ON SAFETY	2
ABOUT BLOOD PRESSURE	6
PRECAUTIONS BEFORE USE	6
FEATURES OF THE PRODUCT	8
PARTS IDENTIFICATION	10
INSERT OR REPLACE BATTERIES	11
MEM. UNIT AND VOICE SYSTEM SETUP	12
UNIT CONVERSION mmHg/kPa DISPLAY	13
WHO BLOOD PRESSURE CLASSIFICATION DISPLAY	14
ATTACHING THE WRIST CUFF	15
HOW TO MEASURE BLOOD PRESSURE	16
CLEAN AND MAINTENANCE	18
SPECIFICATIONS	19
TROUBLESHOOTING	20
CALIBRATION METHOD	29

INTRODUCTION

The Monitor uses the oscillometric method of blood pressure measurement. Measurement Automatic Electronic Blood Pressure Monitor is intended for use by medical professionals or at home to monitor and display diastolic, systolic blood pressure and the pulse rate, with an air wrist cuff buckled around one's wrist according to the instructions in the "ATTACHING THE WRIST CUFF". The expected life of the product is 5 years.

The product complies with the electromagnetic compatibility requirement of IEC 60601-1-2 and safety standards of IEC 60601-1 and performance of IEC 60601-1-2:30 as specified in Regulation (EU)2017/745.

NOTES ON SAFETY

* The warning signs and notes icons shown here are listed for your safe and correct use of the unit, so as to prevent injuries or damages to the device.
* The icons and meanings are as follow.

Examples of signs
The icon indicates prohibitions (what you should not do). Matters involving actual prohibitions are indicated by text or pictures in or near . The left icon refers to "general prohibition".

The icon indicates something that is compulsory (what must always be observed). Matters involving actual compulsory actions are indicated by text or pictures in or near . The left icon refers to "general compulsion".

The icon indicates something can't be disassembled or "Don't disassemble". Matters involving actual compulsory actions are indicated by text or pictures in or near . The left icon refers to "general prohibition".

Type BF Applied part IP Classification: IP20 **Caution**

Please refer to the instructions for use **Consult instruction for use** **The following symbol indicates that the device is MR-unsafe.**

MD Indicates medical device

Indicates a medical device that needs to be protected from moisture.

Contact its local authorities to determine the proper method of disposal of potentially bio-hazardous parts and accessories.

Patient must follow doctor's instruction and should not perform self-judgment and self-treatment by the measuring result. Self-diagnosis of measured results and treatment are dangerous. The device should not be used to judge illness, first aid and continuously monitor measuring.

This device can not be used for Patient transport and surgical care. It can be used in household or fixed places only. Please press "on/off" button to stop work when you feel uncomfortable with the wrist, or if the air is inflating abnormally without stop.

This device should not be used by children under 18 years old or people who cannot express their will, otherwise it will cause harm.

Do not use the unit for purpose other than measuring blood pressure. May cause accident or trouble.

Please do not use mobile phone around the device. Please do not use the device around the magnetic field.

The device is prohibited from being used during movement.

Do not use the equipment in outdoor or shower rooms.

Do not disassemble, repair, or remodel the main unit or the wrist cuff of the blood pressure monitor.

Will cause the unit to function erroneously.

Requests from Manufacturer

Make sure there is no connection tubing kinking before starting measurement to avoid any injury to patient.

For any patient, do not measure more than 3 times continuously, it should be at least above 5 minutes of interval rest between any two measurements, otherwise will cause extravasated blood.

Do not measure your blood pressure over 6 times each day.

Do not apply the cuff over a wound as this can cause further injury.

Do not measure on the wrist which is on the side of a mastectomy, otherwise it could cause injury.

Observe the air pressure value from the LCD display.

When measuring, it could not exceed 280 mmHg, otherwise Please press "on/off" button to stop.

Do not use force to bend the wrist cuff or the air tube.

Do not knock or drop the main unit.

Always use the specified accessories in the manual, the use of other parts not approved by the manufacturer may cause faults or injuries.

For service information, parts list etc., please contact the dealer.

ABOUT BLOOD PRESSURE

1. What is blood pressure?
Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts. Diastolic pressure occurs when the heart expands.

Blood pressure is measured in millimeters of mercury (mmHg). One's natural blood pressure is represented by the fundamental pressure, which is measured first thing in the morning while one is still at rest and before eating.

2. What is hypertension and how is it controlled?
Hypertension, an abnormally high arterial blood pressure, if left unattended, can cause many health problems including stroke and heart attack. Hypertension can be controlled by altering lifestyle, avoiding stress and with medication under a doctor's supervision while you take your medicine. To prevent hypertension or keep it under control:

• Do not smoke
• Reduce salt and fat intake
• Maintain proper weight

• Exercise regularly
• Have regular physical checkups

3. Why measure blood pressure at home?

Blood pressure measured at a clinic or doctor's office may cause apprehension and produce an elevated reading, 25 to 30 mmHg higher than that measured at home. Home measurement reduces the effects of outside influences on blood pressure readings, supplements the doctor's readings and provides a more accurate, complete blood pressure history.

4. WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), and shown in chart below.

5. Blood pressure variations

An individual's blood pressure varies greatly on a daily and seasonal basis. It may vary by 30 to 50 mmHg due to various conditions during the day. In

The PATIENT is an intended OPERATOR.
Not servicing and maintenance while the ME EQUIPMENT is in use.
The user can maintain the product, the maintenance method is described in the maintenance instructions of manual.
Stop using the equipment immediately, if it is in contact with water.

ABOUT BLOOD PRESSURE

1. What is blood pressure?
Blood pressure is the force exerted by blood against the walls of the arteries. Systolic pressure occurs when the heart contracts. Diastolic pressure occurs when the heart expands.

Blood pressure is measured in millimeters of mercury (mmHg). One's natural blood pressure is represented by the fundamental pressure, which is measured first thing in the morning while one is still at rest and before eating.

2. What is hypertension and how is it controlled?
Hypertension, an abnormally high arterial blood pressure, if left unattended, can cause many health problems including stroke and heart attack. Hypertension can be controlled by altering lifestyle, avoiding stress and with medication under a doctor's supervision while you take your medicine. To prevent hypertension or keep it under control:

• Do not smoke
• Reduce salt and fat intake
• Maintain proper weight

• Exercise regularly
• Have regular physical checkups

3. Why measure blood pressure at home?

Blood pressure measured at a clinic or doctor's office may cause apprehension and produce an elevated reading, 25 to 30 mmHg higher than that measured at home. Home measurement reduces the effects of outside influences on blood pressure readings, supplements the doctor's readings and provides a more accurate, complete blood pressure history.

4. WHO blood pressure classification

Standards for assessment of high blood pressure, without regard to age, have been established by the World Health Organization (WHO), and shown in chart below.

5. Blood pressure variations

An individual's blood pressure varies greatly on a daily and seasonal basis. It may vary by 30 to 50 mmHg due to various conditions during the day. In

hypertensive individuals, variations are even more pronounced. Normally, the blood pressure rises while at work or play and falls to its lowest levels during sleep. So, do not be overly concerned by the results of one measurement. Take measurements at the same time every day using the procedure described in this manual, and know your normal blood pressure. Many readings give a more comprehensive blood pressure history. Be sure to note date and time when recording your blood pressure. Consult your doctor to interpret your blood pressure data.

FEATURES OF THE PRODUCT

- Memory can store 60/60 measurements.
- Large and clear LCD display.

PRECAUTIONS BEFORE USE

- If you are taking medication, consult with your doctor to determine the most appropriate time to measure your blood pressure. NEVER change a prescribed medication without first consulting with your doctor.

INSERT OR REPLACE BATTERIES

- Remove the battery cover.
- Insert new batteries into the battery compartment as shown, taking care that the polarities(+) and (-) are correct.
- Close the battery cover. Use only LR03, AAA batteries.

Disposal of empty battery to the authorized collecting party subject to the regulation of each individual territory.

CAUTION
Insert the batteries as shown in the battery compartment. If not, the device will not work.
When (LOW BATTERY mark) blinks in the display, replace all batteries with new ones. Do not mix old and new batteries. It may shorten the battery life, or cause the device to malfunction.
 (LOW BATTERY mark) does not appear when the batteries run out.
Please ensure to distinguish positive polar "+" and negative polar "-" of batteries when replacing batteries.

PARTS IDENTIFICATION

UNIT CONVERSION mmHg/kPa DISPLAY

The goods have mm Hg(mmHg), kPa (kPa) two kinds of blood pressure display units(mmHg factory to express).

WHO BLOOD PRESSURE CLASSIFICATION DISPLAY

Diastolic blood pressure
Reference material: journal of hypertension 1999, vol 17 No.2

- Grade 3 hypertension (severe)
- Grade 2 hypertension (moderate)
- Grade 1 hypertension (mild)
- High-normal
- Normal
- Optimal

ATTACHING THE WRIST CUFF

1. Fastening the wrist cuff
1) Wrap the wrist cuff around your wrist about 1-2cm above your hand as shown in the figure at the right.
2) Fasten the wrist cuff tightly by using the Velcro Strip. Note:
For proper measurements, fasten the wrist cuff tightly and measure on a bare wrist.

2. How to take proper measurements
For best accuracy in blood pressure measurement:
• Sit comfortably at a table. Rest your wrist on the table.
• Relax for about 5 to 10 minutes before measurement.
• Raise your hand so that the wrist cuff is at the same level as your heart.
• Remain still and keep quiet during measurement.
• Do not measure right after physical exercise or a bath.
• Rest twenty or thirty minutes before taking measurement every day.
• Measure your blood pressure at about the same time every day.

HOW TO MEASURE BLOOD PRESSURE

- Set up the wrist cuff to your wrist as previous section of "ATTACHING THE WRIST CUFF".
- Press or button, the display shows "0". Inflation icon flashes on display when it begins to inflate.
- During the inflation, please do not move
- Automatic deflation after measurement and display the blood pressure, heart rate, and blood pressure indicator, voice report (monitor with talking system only)
- Mistake, display "E", please refer to the instruction about troubleshooting.
- If cuff remove during the measurement, please fix the cuff and test again.
- Stop the measurement in emergency, please press or to cut off the power.
- Automatic deflation at the same time.
- Power Off, Press or to cut off the power.

MEMORY STORAGE

- under power off, hold or for 3 seconds, get into relative memory mode, it reads average of latest 3 measurements, press for the buttons (UP), "SET" button for the memory (DOWN)
* voice report (monitor with talking system only)
* After measurement, press "SET"
- Press or to cut off the power.
* If you forget to turn it off, will automatically power off in one minute.
* More than 60 measurements, it will deletes the earliest one automatically.
* The stored memory will not lose after replacing battery.
- Delete memory (all the memory will be deleted) Showing memory, hold "SET" for 3 seconds, display "00" means all the memories have been deleted.

CLEAN AND MAINTENANCE

- Keep this device in the case provided with the device when you do not use it.
- Do not fold the arm cuff too tightly.

SPECIFICATIONS

Measuring Method	Oscillometric Measurement
Indication	Digital LCD display
Measuring Range:	Pressure: (30-280)mmHg Pulse: (40-199)Beat/min
Accuracy:	Static Pressure: ± 3mmHg Pulse: ± 5%
Memory:	90 Memories
Power supply:	2x1.5V Batteries(LR3 or AAA) use alkaline battery, measure above 200 times.
Operating condition:	+5°C~+40°C, 30%RH~80%RH Atmospheric pressure: 86kPa~106kPa
Storage condition:	-20°C~+55°C, 10%RH~93%RH Atmospheric pressure: 50kPa~106kPa
Dimensions:	Approx: 78(W)X78(H)X34(D)mm
Weight:	Approx: 125g, excluding batteries
Classification	Type BF
Wrist circumference	(13.5-19.5)cm

* Specifications may be changed without notice in the event of improvement being made.

TROUBLESHOOTING

If you have trouble in using the unit please check the following points first.

ERROR DISPLAY	POSSIBLE CAUSE	HOW TO CORRECT
Nothing is displayed	No battery installation	Insert batteries
When you push the POWER button or BATTERY icon flash	Battery worn out The polarities of batteries placed wrongly	Replace new batteries Insert battery in the correct polarities

Guidance and manufacturer's declaration – electromagnetic immunity

The Model PG-800AD Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800AD Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	+8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60Hz	30 A/m, 50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_i is the a.c. mains voltage prior to application of the test level

Guidance and manufacturer's declaration – electromagnetic immunity

The Model PG-800AD Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800AD Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 6 Vrms 80 MHz to 150 kHz 150 kHz outside ISM bands	N/A	Portable and mobile RF communications equipment should be used no closer to any part of the Model PG-800AD Series Electronic Blood Pressure Monitor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \frac{3.5}{f} \sqrt{P}$

Guidance and manufacturer's declaration – electromagnetic immunity

The Model PG-800AD Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800AD Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	+8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60Hz	30 A/m, 50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_i is the a.c. mains voltage prior to application of the test level

1. Degree of protection against electric shock: INTERNALLY POWERED EQUIPMENT.
2. Degree of protection against electric shock: TYPE BF APPLIED PART.
3. Mode of operation: CONTINUOUS OPERATION.
4. Equipment not suitable for category AP&APG equipment use in presence.

STATEMENT
The system might not meet its performance specifications if stored or used outside the environment and humidities as mentioned below.
Operating conditions: +5°C~+40°C, 30%RH~80%RH 86kPa~106kPa
Storage conditions: -20°C~+55°C, 10%RH~93%RH

Guidance and manufacturer's declaration – electromagnetic immunity

The Model PG-800AD Series Electronic Blood Pressure Monitor is intended for use in the electromagnetic environment specified below. The customer or the user of the Model PG-800AD Series Electronic Blood Pressure Monitor should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	+8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m, 50/60Hz	30 A/m, 50/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE U_i is the a.c. mains voltage prior to application of the test level

1. Degree of protection against electric shock: INTERNALLY POWERED EQUIPMENT.
2. Degree of protection against electric shock: TYPE BF APPLIED PART.
3. Mode of operation: CONTINUOUS OPERATION.
4. Equipment not suitable for category AP&APG equipment use in presence.

STATEMENT
The system might not meet its performance specifications if stored or used outside the environment and humidities as mentioned below.
Operating conditions: +5°C~+40°C, 30%RH~80%RH 86kPa~106kPa
Storage conditions: -20°C~+55°C, 10%RH~93%RH

Calibration Method

- Press and hold the "ON/OFF, MEM" button at the same time, load the battery, enter the static air pressure calibration mode after the LCD screen is fully displayed, and then release the button.
- Press ON/OFF to close the internal air valve.
- Connect the external standard barometric interface and the digital barometer interface to the cuff interface.

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a The ISM (industrial, scientific and medical) bands between 0.15 MHz and 80 MHz are 6.765 MHz to 6.975 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.660 MHz to 40.70 MHz. The amateur radio bands between 0.15 MHz and 80 MHz are: 1.8 MHz to 2.0 MHz, 3.5 MHz to 4.0 MHz, 5.3 MHz to 5.4 MHz, 7 MHz to 7.3 MHz, 10.1 MHz, 10.15 MHz, 14.1 MHz to 14.2 MHz, 18.07 MHz to 18.17 MHz, 21.0 MHz to 21.4 MHz, 24.89 MHz to 24.99 MHz, 28.0 MHz to 29.7 MHz and 50.0 MHz to 54.0 MHz.

b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.7 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formula used in calculating the recommended separation distance for transmitters in these frequency ranges.

For transmitters rated at a maximum output power not listed above the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

CALIBRATION METHOD

- Press and hold the "ON/OFF, MEM" button at the same time, load the battery, enter the static air pressure calibration mode after the LCD screen is fully displayed, and then release the button.
- Press ON/OFF to close the internal air valve.
- Connect the external standard barometric interface and the digital barometer interface to the cuff interface.

4. External input 50mmHg and 200mmHg standard static air pressure, and observe the air pressure value displayed at the position of the LCD systolic pressure (SYS) and the value of the digital pressure gauge should be in the range of +/-3mmHg.

Caution
1. ME devices can be used in exposed environments, including electromagnetic interference environment to ensure basic safety and basic performance unchanged. In the event of any serious event related to this product, such as a serious adverse event, significant alteration of the product resulting in change of intended use, etc., it will be reported to the manufacturer and the competent authorities of the user and/or the member states where the patient is located.

Notes:
Essential performance: Limits of error of the manometer, ±3mmHg. Reproducibility of the blood pressure determination, ±3mmHg.
Clinical benefits: Accurate measurement with SBP and DBP, clinical performance meets the requirements of ISO 81060-2:2018.